

Instructions μ-Slide 8 Well high



The ibidi product family is comprised of a variety of μ -Slides, μ -Dishes, and μ -Plates which have all been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength. The μ -Slide 8 Well high is an 8 well chamber slide in which cells can be cultivated and, subsequently, investigated with microscopical methods. This open μ -Slide (chambered coverslip) is intended for cell culture, immunofluorescence, live cell imaging, and high-end microscopy.

Overview

This document is applicable to the following product numbers:

| Cat. No. | Product Name |
|----------|---|
| 80806 | μ-Slide 8 Well ^{high} ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized, individually packed |
| 80806-90 | μ -Slide 8 Well high ibiTreat, Bulk Pack: #1.5 polymer coverslip, tissue culture treated, sterilized, individually packed |
| 80806-96 | μ -Slide 8 Well high ibiTreat, Bulk Pack: #1.5 polymer coverslip, tissue culture treated, sterilized, 8 per tray, 12 trays |
| 80809 | μ-Slide 8 Well ^{high} Collagen I: #1.5 polymer coverslip, sterilized, individually packed |
| 80802 | μ-Slide 8 Well ^{high} Collagen IV: #1.5 polymer coverslip, sterilized, individually packed |
| 80804 | μ-Slide 8 Well ^{high} Poly-L-Lysine: #1.5 polymer coverslip, sterilized, individually packed |
| 80801 | μ -Slide 8 Well $^{\rm high}$ Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized, individually packed |

Material

ibidi μ -Slides, μ -Dishes, and μ -Plates are made of a polymer that has the highest optical quality. The polymer coverslip on the bottom exhibits extremely low birefringence and autofluorescence, similar to that of glass. Also, it is not possible to detach the bottom from the upper part. The μ -Slides, μ -Dishes, and μ -Plates are intended for one-time use and are not autoclavable, since they are only temperature-stable up to 80° C/175°F. Please note that gas exchange between the medium and the incubator's atmosphere occurs partially through the polymer coverslip, which should not be covered.

| Optical Properties ibidi Polymer Coverslip | | |
|--|-------------------|--|
| Refractive index n _D (589 nm) | 1.52 | |
| Abbe number | 56 | |
| Thickness | No. 1.5 (180 μm) | |
| Material | Polymer coverslip | |

Please note! The ibidi Polymer Coverslip is compatible with certain types of immersion oil only. A list of suitable oils can be found on page 3.

Shipping and Storage

The μ -Slides, μ -Dishes and μ -Plates are sterilized and welded in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is listed in the following table.


| Conditions | | | |
|-------------------------|--------------|--------------|--|
| Shipping conditions | Ambient | | |
| Storage conditions | RT (15–25°C) | RT (15–25°C) | |
| Shel | f Life | | |
| | | | |
| ibiTreat, Uncoated | 36 months | | |
| Collagen I, Collagen IV | 18 months | | |
| Poly-L-Lysine | 18 months | | |

Geometry of the μ -Slide 8 Well high

The $\mu\text{-Slide}$ 8 Well high provides a standard slide format according to ISO 8037/1.



Instructions μ-Slide 8 Well high



| Geometry | | |
|---|--------------------------|--|
| Outer dimensions in mm $(w \times l)$ | 25.5 × 75.5 | |
| Number of wells | 8 | |
| Dimensions of wells in mm $(w \times l \times h)$ | $9.4\times10.7\times9.3$ | |
| Volume per well | 300 µl | |
| Height with/without lid | 10.8/9.5 mm | |
| Growth area per well | 1.0 cm^2 | |
| Coating area per well | $2.2\mathrm{cm}^2$ | |
| Bottom | ibidi Polymer Coverslip | |

Surface

The tissue culture-treated ibiTreat surface is a physical surface modification and optimized for adhesion of most cell types. The uncoated surface is a very hydrophobic surface and allows no direct cell growth. It is suitable for specific coatings or suspension cells.

If you like to establish a particular coating for your demands we recommend testing your coating procedure on uncoated and ibiTreat surfaces, since some proteins and biomolecules adhere differently to hydrophobic or hydrophilic polymer surfaces.

The μ-Slide 8 Well high is also available with a Collagen Type I, Collagen Type IV or a Poly-L-Lysine coated surface. For the coating, only high quality proteins are used: Collagen Type I: ibidi #50203, Collagen Type IV: Corning #356233, Poly-L-Lysine: Sigma #P4832.

Coating

Detailed information about coatings is provided in Application Note 08: Coating protocols for ibidi labware products.

In short, specific coatings are possible following this protocol:

1. Prepare your coating solution according to the manufacturer's specifications or reference.

- 2. Apply 300 µl and leave at room temperature for at least 30 minutes.
- 3. Aspirate the solution and wash with the recommended protein dilution buffer.
- 4. The μ -Slide 8 Well ^{high} is ready to be used. Optionally let dry at room temperature. Attention, some coating proteins might degenerate when drying!

Seeding Cells

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, application of a $5-11 \times 10^4$ cells/ml suspension should result in a confluent layer within 2-3 days.
- Apply 300 µl cell suspension into each well. Avoid shaking as this will result in inhomogeneous distribution of the cells.
- Cover the slide with the supplied lid. Incubate at 37°C and 5 % CO₂ as usual.

Undemanding cells can be left in their seeding medium for up to three days and grow to confluence there. However, best results might be achieved when the medium is changed every 1–2 days. Carefully aspirate the old medium and replace it by 300 µl fresh medium per well.

Tip:

As you may know from 96 well plates, the bent meniscus at the air–liquid interphase in small open wells destroys the phase contrast effect of your microscope image. To avoid this problem, we recommend using our channel Slides such as the $\mu\text{-Slides}$ I Luer and $\mu\text{-Slide}$ VI $^{0.4}$ or a Ph+ Slide.



Instructions μ-Slide 8 Well high

Chemical Compatibility

The following table provides some basic information on the chemical and solvent compatibility of the μ -Slide 8 Well ^{high}. For a full list of compatible solvents and more information on chemical compatibility, please visit the FAQ section on ibidi.com.

| Chemical / Solvent | Compatibility |
|--------------------|-------------------------------------|
| Methanol | yes |
| Ethanol | yes |
| Formaldehyde | yes |
| Acetone | yes, without lid |
| Mineral oil | no |
| Silicone oil | yes |
| Immersion oil | See Immersion Oil on page 3. |

Microscopy

To analyze your cells, no special preparations are necessary. Cells can be directly observed live or fixed, preferably on an inverted microscope. The bottom cannot be removed. For optimal results in fluorescence microscopy and storage of fixed and stained samples, ibidi provides mounting media (50001 and 50011) optimized for $\mu\text{-Dishes}, \mu\text{-Slides},$ and $\mu\text{-Plates}.$

Immersion Oil

When using oil immersion objectives with the ibidi Polymer Coverslip, use only the immersion oils specified in the table below. The use of any non-recommended oil could damage the ibidi Polymer Coverslip. The resulting leakage may harm objectives and microscope components. All immersion oils that are not listed in the table below should be considered as non-compatible.

| Company | Product | Ordering No. | Lot Number | Test Date |
|-----------|-----------------------------|---------------|------------|-----------|
| ibidi | ibidi Immersion Oil | 50101 | 16-12-27 | 01/2017 |
| Cargille | Type A | 16482 | 100592 | 01/2017 |
| Cargille | Type HF | 16245 | 92192 | 01/2017 |
| Carl Roth | Immersion oil | X899.1 | 414220338 | 01/2017 |
| Leica | Immersion Liquid | 11513859 | n.a. | 03/2011 |
| Nikon | Immersion Oil F2 30cc | MXA22192 | n.a. | 01/2020 |
| Nikon | Silicone Immersion Oil 30cc | MXA22179 | 20191101 | 01/2020 |
| Olympus | Silicone Immersion Oil | SIL300CS-30CC | N4190800 | 01/2017 |
| Zeiss | Immersol 518 F | 444960 | 160706 | 01/2017 |
| Zeiss | Immersol W 2010 | 444969 | 101122 | 04/2012 |



μ-Slide 8 Well high

Instructions

Ordering Information

μ-Slide 8 Well high



| Cat. No. | Description |
|----------|--|
| 80806 | μ-Slide 8 Well ^{high} ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized, individually packed |
| 80806-90 | μ-Slide 8 Well ^{high} ibiTreat, Bulk Pack : #1.5 polymer coverslip, tissue culture treated, sterilized, individually packed |
| 80806-96 | μ-Slide 8 Well ^{high} ibiTreat, Bulk Pack : #1.5 polymer coverslip, tissue culture treated, sterilized, 8 per tray, 12 trays |
| 80809 | μ-Slide 8 Well ^{high} Collagen I: #1.5 polymer coverslip, sterilized, individually packed |
| 80802 | μ-Slide 8 Well ^{high} Collagen IV: #1.5 polymer coverslip, sterilized, individually packed |
| 80804 | μ-Slide 8 Well ^{high} Poly-L-Lysine : #1.5 polymer coverslip, sterilized, individually packed |
| 80801 | μ-Slide 8 Well ^{high} Uncoated : #1.5 polymer coverslip, hydrophobic, sterilized, individually packed |
| 80800 | μ-Slide 8 Well ^{high} Bioinert : #1.5 polymer coverslip, surface passivation with Bioinert, sterilized, individually packed |
| 80807 | μ -Slide 8 Well ^{high} Glass Bottom: #1.5H (170 μ m ±5 μ m) D 263 M Schott glass, sterilized, individually packed |
| 80807-90 | μ-Slide 8 Well ^{high} Glass Bottom, Bulk Pack : #1.5H (170 μm ±5 μm) D 263 M Schott glass, sterilized, individually packed |
| 80807-96 | μ-Slide 8 Well ^{high} Glass Bottom, Bulk Pack : #1.5H (170 μm ±5 μm) D 263 M Schott glass, sterilized, 8 per tray, 12 trays |

μ-Slide 8 Well high Grid-500



| Cat. No. | Description |
|------------|--|
| 80806-G500 | μ-Slide 8 Well ^{high} ibiTreat Grid-500 : #1.5 polymer coverslip, tissue culture treated, grid repeat distance 500 μm, sterilized, individually packed |

For research use only!

Further information can be found at ibidi.com. For questions and suggestions please contact us by e-mail *info@ibidi.de* or by telephone +49 (0)89/520 4617 0.

 $\ensuremath{\mathbb{C}}$ ibidi GmbH, Lochhamer Schlag 11, 82166 Gräfelfing, Germany.