

The ibidi product family is comprised of a variety of μ -Slides and μ -Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength.

The μ -Slide VI ^{0.1} is designed for flow assays in a minimal volume. It can be connected to a pump and enables you to observe cells under flow conditions. The small dimensions offer you the possibility to work with a minimum of cells (e.g. mouse model). The convenient six channel format is ideal for the application of standard protocols for e.g. treatment, staining, and microscopy of living or fixed cells.

Material

ibidi μ -Slides, μ -Dishes, and μ -Plates are made of a polymer that has the highest optical quality. The polymer coverslip on the bottom exhibits extremely low birefringence and autofluorescence, similar to that of glass. Also, it is not possible to detach the bottom from the upper part. The μ -Slides, μ -Dishes, and μ -Plates are intended for one-time use and are not autoclavable, since they are only temperature-stable up to 80° C/175°F. Please note that gas exchange between the medium and the incubator's atmosphere occurs partially through the polymer coverslip, which should not be covered.

Optical Properties ibidi Polymer Coverslip		
Refractive index n _D (589 nm)	1.52	
Abbe number	56	
Thickness	No. 1.5 (180 μm)	
Material	Polymer coverslip	

Please note! The ibidi Polymer Coverslip is compatible with certain types of immersion oil only. A list of suitable oils can be found on page 3.

Shipping and Storage

The μ -Slides, μ -Dishes and μ -Plates are sterilized and welded in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is listed in the following table.


Conditions		
Shipping conditions	Ambient	
Storage conditions	RT (15–25°C)	

Shelf Life		
ibiTreat, Uncoated	36 months	
Collagen IV	18 months	

Geometry

The μ -Slide VI ^{0.1} provides a standard slide format according to ISO 8037/1. The lateral adapter to adapter distance of 9 mm (like 96 well plates) allows using multichannel pipettes.

Geometry of the μ -Slide VI $^{0.1}$		
Outer dimensions	25.5 mm x 75.5 mm	
Adapters	Female Luer	
Number of channels	6	
Channel volume	1.7 µl	
Channel height	0.1 mm	
Channel length	17 mm	
Channel width	1.0 mm	
Volume per reservoir	60 µl	
Growth area	0.17 cm ² per channel	
Coating area using 1.7 µl	0.34 cm ² per channel	
Bottom	No. 1.5 ibidi Polymer Cov- erslip	





Surface

The tissue culture-treated ibiTreat surface is a physical surface modification and optimized for adhesion of most cell types. The uncoated surface is a very hydrophobic surface and allows no direct cell growth. It is suitable for specific coatings or suspension cells.

If you like to establish a particular coating for your demands we recommend testing your coating procedure on uncoated and ibiTreat surfaces, since some proteins and biomolecules adhere differently to hydrophobic or hydrophilic polymer surfaces.

Coating


Detailed information about coatings is provided in Application Note 08: Coating protocols for ibidi labware products.

In short, specific coatings are possible following this protocol:

- 1. Prepare your coating solution according to the manufacturer's specifications or reference.
- 2. Apply 1.7 µl and leave at room temperature for at least 30 minutes.
- 3. Aspirate the solution and wash with the recommended protein dilution buffer.
- 4. The μ-Slide VI ^{0.1} is ready to be used. Optionally let dry at room temperature. Attention, some coating proteins might degenerate when drying!

Filling and Handling

Filling the very small channels of μ -Slide VI ^{0.1}, especially hydrophobic uncoated, by a normal pipette might be challenging. Use a 20 μ l pipet and directly inject into the channel opening. Optionally, use a small volume syringe with a Luer tip for convenient filling the channels with coating solution or cell suspension.



Seeding Cells

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, application of a $12-28 \times 10^5$ cells/ml suspension should result in a confluent layer within 2–3 days.
- Apply 1.7 μl cell suspension into the channel of the μ-Slide. Quick dispensing helps to avoid trapped air bubbles.
- Cover reservoirs with the supplied lid. Incubate at 37° C and 5% CO₂ as usual.
- Await cell attachment in order not to flush out the cells. Afterwards fill each reservoir with 60 μ l cell free medium.
- Connect the µ-Slide to the pump and conduct your perfusion experiment.

Tip:

The day before seeding the cells we recommend placing the cell medium and the μ -Slide into the incubator for equilibration. This will prevent the liquid inside the channel from emerging air bubbles over the incubation time.

Important!

After coating the µ-Slide with a coating that must not be dried, seed cells without emptying the channel: First, aspirate all remaining liquid from both reservoirs. Do not empty the channel. Then, fill 6 µl of cell suspension onto one channel opening inside a Luer reservoir. After that, slowly remove 6 µl from the opposite Luer reservoir. Make sure to avoid trapped air bubbles.



Exchanging Medium

Aspirate both reservoirs and slowly fill 60 μ l of fresh medium into each of the reservoirs. The equilibration of the liquid levels may take some minutes because of the small width of the channel.

Important!

Please note that the μ -Slide VI ^{0.1} is not for use in static culture due to the small channel volume. Cultivation without perfusion is only possible by a medium exchange every few hours or by using a rocker which constantly generates a slight medium flow between the two reservoirs.

Microscopy

To analyze your cells, no special preparations are necessary. Cells can be directly observed live or fixed, preferably on an inverted microscope. The bottom cannot be removed. For optimal results in fluorescence microscopy and storage of fixed and stained samples, ibidi provides a mounting medium (50001) optimized for μ -Dishes, μ -Slides, and μ -Plates.

Chemical Compatibility

The following table provides some basic information on the chemical and solvent compatibility of the μ -Slide VI ^{0.1}. For a full list of compatible solvents and more information on chemical compatibility, please visit the FAQ section on ibidi.com.

Compatibility
yes
yes
yes
yes, without lid
no
yes
See Immersion Oil on page 3.

Immersion Oil

When using oil immersion objectives with the ibidi Polymer Coverslip, use only the immersion oils specified in the table below. The use of any non-recommended oil could damage the ibidi Polymer Coverslip. The resulting leakage may harm objectives and microscope components. All immersion oils that are not listed in the table below should be considered as non-compatible.

Company	Product	Ordering No.	Lot Number	Test Date
ibidi	ibidi Immersion Oil	50101	16-12-27	01/2017
Cargille	Type A	16482	100592	01/2017
Cargille	Type HF	16245	92192	01/2017
Carl Roth	Immersion oil	X899.1	414220338	01/2017
Leica	Immersion Liquid	11513859	n.a.	03/2011
Nikon	Immersion Oil F2 30cc	MXA22192	n.a.	01/2020
Nikon	Silicone Immersion Oil 30cc	MXA22179	20191101	01/2020
Olympus	Silicone Immersion Oil	SIL300CS-30CC	N4190800	01/2017
Zeiss	Immersol 518 F	444960	160706	01/2017
Zeiss	Immersol W 2010	444969	101122	04/2012




Instructions

Ordering Information


The $\mu\text{-Slide VI}$ family is available in different surfaces and bottom characteristics. $\mu\text{-Slide VI}^{0.4}$

	Cat. No.	Description
	80606	μ-Slide VI ^{0.4} ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
60000	80606-90	μ -Slide VI ^{0.4} ibiTreat, Bulk Pack: #1.5 polymer coverslip, tissue culture treated, sterilized
13 00000	81602	μ-Slide VI ^{0.4} Collagen IV: #1.5 polymer coverslip, sterilized
	81604	μ-Slide VI ^{0.4} Poly-L-Lysine : #1.5 polymer coverslip, sterilized
	81601	μ -Slide VI ^{0.4} Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized


$\mu\text{-Slide VI}^{\,0.5}$ Glass Bottom



µ-Slide VI ^{0.1}



Tube Adapter Set



For research use only!

Further information can be found at www.ibidi.com. For questions and suggestions please contact us by e-mail *info@ibidi.de* or by telephone +49 (0)89/520 4617 0.

© ibidi GmbH, Lochhamer Schlag 11, 82166 Gräfelfing, Germany.