

ibidi Application Guide

3D Cell Culture

3D Cell Culture	2
Matrices for 3D Cell Culture	3
3D Cell Culture Assays	4
Spheroid and Organoid Culture	4
Single Cells in a 3D Matrix	6
Chemotaxis and Migration Assays in 3D	7
3D Cell Culture Under Flow	8
Experimental Examples	9
Single Cells in a 3D Matrix	9
Spheroid and Organoid Culture	9
Chemotaxis and Migration Assays in 3D	11

Selected Publications

P. Xu et al. Modulation of Intestinal Epithelial Permeability by Plasma from Patients with Crohn's Disease in a Three-dimensional Cell Culture Model. Scientific Reports, 2019, 10.1038/s41598-018-38322-8

[Read article](#)

E. Hoque Apu, S.U. Akram, J. Rissanen, H. Wan and T. Salo. Desmoglein 3 – Influence on oral carcinoma cell migration and invasion. Experimental Cell Research, 2018, 10.1016/j.yexcr.2018.06.037

[Read article](#)

M. Dietrich et al. Guiding 3D cell migration in deformed synthetic hydrogel microstructures. Soft Matter, 2018, 10.1039/C8SM00018B

[Read article](#)

H. Grobe, A. Wüstenhagen, C. Baarlink, R. Grosse and K. Grikscheit. A Rac1-FMNL2 signaling module affects cell-cell contact formation independent of Cdc42 and membrane protrusions. PloS one, 2018, 10.1371/journal.pone.0194716

[Read article](#)


3D Cell Culture

The majority of cells in living tissue grow in a three-dimensional microenvironment, where they communicate and interact with each other and their surroundings. Animal cells are embedded in the extracellular matrix (ECM), which is composed of proteoglycans and fibrous proteins (mainly collagen, elastin, and fibronectin). This complex, dynamic, and tissue-specific 3D structure provides physical scaffolding for the cells and initiates cues that influence cell differentiation and behavior.


When cultured in a traditional, two-dimensional *in vitro* environment, the cells are attached to a flat surface (e.g., a monolayer in a standard cell culture dish) and can only grow and migrate on the substrate. In a 3D *in vitro* setup, the cells are grown in suspension on a non-adhesive surface, or they can be embedded in or on a 3D matrix (e.g., Matrigel® or collagen I) that mimics the ECM and allows them to grow in all three directions.

Cells behave differently inside a 3D gel matrix compared to a 2D environment. In many cases, a 3D environment reflects the *in vivo* situation more accurately. This should be considered when analyzing cell behavior, differentiation, response to drug treatment, and gene and protein expression.

Not surprisingly, many cell culture approaches have been adapted to a 3D environment. This includes drug screenings that use [spheroids and organoids](#), which are indispensable nowadays as tumor models.



Confocal laser scanning microscopy projection of an HT-1080 LifeAct spheroid. The colors indicate the distance from the surface. Warm colors = close to the surface, cold colors = distant from the surface.



3D cell culture applications in comparison to the traditional 2D monolayer.

Edmondson R, Broglie JJ, Adcock AF, Yang L (2014) Three-dimensional cell culture systems and their applications in drug discovery and cell-based biosensors. *Assay Drug Dev Technol* 12(4):207–18. 10.1089/adt.2014.573.

[Read article](#)

Frantz C, Stewart KM, Weaver VM (2010) The extracellular matrix at a glance. *J Cell Sci* 123(Pt 24):4195–200. 10.1242/jcs.023820.

[Read article](#)

Matrices for 3D Cell Culture

A variety of matrices are currently being used for culturing cells in a 3D environment, including a wide range of natural proteins to synthetic scaffolds. The choice of a suitable matrix strongly depends on the cell type being used and specific experimental questions. Hydrogels, such as collagen and fibrin, are water-swollen polymer networks that are commonly used. They can substitute the natural matrix in 3D cell culture experiments because they mimic several key features of the native extracellular matrix (ECM).


Collagen I is the main component of connective tissue and is abundant in the mammalian body. The fibrous protein consists of three α -chains that combine to create a rope-like triple helix, thus providing tensile strength to the ECM. The triple helices aggregate and form fibrils in a self-organized manner. *In vivo*, the fibrils polymerize into fibers to form tissue such as tendons or dermis.

Collagen (mostly type I) is widely used in 3D cell culture for modeling the extracellular matrix. To embed the cells in the matrix, they are mixed with the liquid gel and pipetted into the vessel. Raising the pH and temperature leads to self-assembly of the collagen fibrils, which results in gelation and the encapsulation of the cells.

Matrigel® is a collagen- and laminin-containing hydrogel used in many 3D cell culture approaches (e.g., organoid culture). It is derived from Engelbreth–Holm–Swarm (EHS) mouse sarcoma tumors and consists mainly of laminin, type IV collagen, and entactin.

Detailed comparison of hydrogels, including their applications, advantages, and disadvantages:

Caliri SR, Burdick JA (2016) A practical guide to hydrogels for cell culture. Nat Methods 13(5):405–414. 10.1038/nmeth.3839
[Read article](#)



The structure of collagen.

ibidi Solution

The ibidi **Collagen Type I, Rat Tail** is a non-pepsinized, native collagen for modeling ECM in gel matrices. Its fast polymerization facilitates optimal cell distribution in 3D gels.

Read our Application Note about how to prepare a 3D gel using the ibidi Collagen Type I, Rat Tail here: [AN 26: Collagen I Gel for 3D Cell Culture](#).



3D Cell Culture Assays

Spheroid and Organoid Culture

Spheroids are cells that adhere to each other under three-dimensional, non-adherent culture conditions. They lack stem cells, which means that they consist of fully differentiated cells. They can be generated by placing them into a scaffold-free suspension using the hanging drop or forced floating method, for example.

Spheroids are not capable of self-renewal and further differentiation. Tumor cell spheroids are an exception, because due to the unlimited proliferation capacities of the tumor cells, they are able to divide and renew. Therefore, spheroids are a useful model for examining tumor cell behavior, such as large-scale drug screenings.

Read here to see a detailed protocol for spheroid generation in the μ -Plate Angiogenesis: [AN 32: Generation of Spheroids \(PDF\)](#)

Organoids are cultured “mini organs”. They can be generated from adult stem cells (ASCs) or pluripotent stem cells (PSCs). When cultured in a three-dimensional matrix/scaffold (e.g., Matrigel® or collagen), these cells differentiate into organ-specific cell types that build small functional organs.

The first generation of intestinal organoids, created from an Lgr5+ stem cell by Sato et al., initiated many protocols for organoid generation from different organs, such as intestine, liver, brain, prostate, kidney, pancreas, lung, and thyroid. Importantly, they can be edited using technologies such as CRISPR, making them a powerful tool for studies about personal therapy, organogenesis, and drug screening.

Sato T, et al. (2009) Single Lgr5 stem cells build crypt-villus structures in vitro without a mesenchymal niche. *Nature* 459(7244):262–265. 10.1038/nature07935.


[Read article](#)

Drost J, Clevers H (2018) Organoids in cancer research. *Nat Rev Cancer* 18:407–418. 10.1038/s41568-018-0007-6.


[Read article](#)

Tuveson D, Clevers H (2019) Cancer modeling meets human organoid technology. *Science* 364(6444):952–955. 10.1126/science.aaw6985.


[Read article](#)



NIH-3T3 cells forming defined spheroids on the [ibidi \$\mu\$ -Pattern](#). Cells were seeded on 200 μ m adhesion spots in a μ -Slide VI^{0.4} and kept under flow (3 dyn/cm²) for 14 days.




Spheroids are cell aggregates, which are often generated from cancer cells.




Organoids are cultured miniature versions of organs, which are derived from stem cells.

ibidi Solutions


The [μ-Slide Spheroid Perfusion](#) is a specialized flow chamber for long-term spheroid culture. Each of the 3 x 7 wells forms its own niche, in which the specimen is cultured. The application of perfusion through the channel on top of the wells ensures optimal nutrition and oxygen diffusion throughout the experiment, without exposing the specimen to significant shear forces.




The [μ-Slides With Multi-Cell μ-Pattern](#) enable spatially defined cell adhesion for spheroid and organoid generation, long-term culture, and high resolution imaging. Defined adhesion spots are able to catch all adherent single cells from a cell suspension. The surrounding Bioinert surface is fully non-cell-attachable. This forces all cells to aggregate to each other at the adhesion spots, thus forming spheroids in a defined and controllable way.




[Bioinert](#) is a stable, biologically inert surface for long-term culture and high-resolution microscopy of spheroids, organoids, and suspension cells on a non-adherent surface without any cell or biomolecule adhesion. It is currently available as the [μ-Dish^{35 mm, high} Bioinert](#), the [μ-Slide 8 Well^{high} Bioinert](#), the [μ-Slide 4 Well Bioinert](#), and the [μ-Slide VI^{9.4} Bioinert](#).




In the [μ-Slide III 3D Perfusion](#), spheroids or organoids can be cultivated in or on a gel layer or embedded in a 3D matrix. The special channel geometry allows for superfusion with a low flow rate (e.g., when utilizing the [ibidi Pump System](#)). This setup makes long-term cultivation possible for up to several weeks. Additionally, the thin coverslip bottom allows for high-resolution imaging.



The [μ-Slide Angiogenesis](#) or the [μ-Plate Angiogenesis 96 Well](#) are easy, cost-effective solutions for the 3D cultivation and microscopy of spheroids and organoids on, or in, gel matrices. The gel layer is directly connected to the medium reservoir above, which enables fast and easy medium exchange by diffusion. For special applications, the [μ-Slide Angiogenesis Glass Bottom](#), with a No. 1.5H glass bottom, is also available.




The ibidi [Collagen Type I, Rat Tail](#) is a non-pepsinized, native collagen for modeling ECM in gel matrices. Its fast polymerization facilitates optimal cell distribution in 3D gels. Read our Application Note about how to prepare a 3D gel using the ibidi Collagen Type I, Rat Tail here: [AN 26: Collagen I Gel for 3D Cell Culture](#).



Single Cells in a 3D Matrix

In many cases, a 3D environment more closely resembles an *in vivo* situation than a 2D cell culture. Single cells can be cultured and imaged in a 3D gel in order to analyze diverse biological questions, such as cell deformation, migration, tube formation, or ECM degradation. In addition to cultures with only one cell type, the invasion behavior of two different cell types (e.g., cancer cells and fibroblasts) can also be investigated by co-culturing them in the same vessel.

In order to isolate cells from the gel matrix, the matrix can be degraded enzymatically (e.g., collagen by collagenase). After this, the cells can be either expanded in a new gel matrix or further processed for DNA, RNA, or protein isolation.



LifeAct-expressing HT-1080 cells (green) in a [Collagen Type I, Rat Tail](#) layer in the [μ-Slide Chemotaxis](#).

ibidi Solutions

The ibidi [Collagen Type I, Rat Tail](#) is a non-pepsinized, native collagen for modeling biological ECM in gel matrices. Its fast polymerization facilitates optimal cell distribution in 3D gels. Read our Application Note about how to prepare a 3D gel using the ibidi Collagen Type I, Rat Tail here: [AN 26: Collagen I Gel for 3D Cell Culture](#).

In the [μ-Slide III 3D Perfusion](#), single cells are embedded in a 3D matrix. The special channel geometry allows for superfusion with a low flow rate (e.g., when utilizing the [ibidi Pump System](#)). Unlike in static cultures, the superfusion ensures optimal oxygen and nutrient supply. This setup makes long-term cultivation possible for up to several weeks. Additionally, the thin coverslip bottom allows for high-resolution imaging.

The [μ-Slide Angiogenesis](#) or the [μ-Plate Angiogenesis 96 Well](#) allow for easy, cost-effective cultivation and microscopy of single cells and co-cultures on, or in, 3D gels. The gel layer is directly connected to the medium reservoir above, enabling fast and easy medium exchange by diffusion. For special applications, the [μ-Slide Angiogenesis Glass Bottom](#), with a No. 1.5H glass bottom, is also available.

The [μ-Slide I Luer 3D](#) is designed for culturing cells on, or in, a 3D gel matrix with defined flow. Each of the three wells can be filled with a gel, in which cells can be embedded. For defined flow application, the channel on top can be connected to a pump (e.g., to the [ibidi Pump System](#)) to ensure optimal oxygen and nutrient supply.

The [μ-Slide Chemotaxis](#) and the [sticky-Slide Chemotaxis](#) are ideal for analyzing single cell migration in 2D and 3D. Chemotactic gradients can be easily established in water-based 3D gels, such as [Collagen I](#) gels and Matrigel®, because the gel structure does not hinder the formation of a soluble gradient by diffusion.

The [μ-Slide Membrane ibiPore Flow](#) is specialized for transmigration and transport studies under both static and flow conditions. It is ideally suited for cell polarity assays, where factors inside a 3D gel matrix lead to the polarization of a cell monolayer. Also, the rolling, adhesion, and transmigration of leukocytes towards chemoattractant-producing cancer cells in a 3D matrix can be observed in high-resolution.

Most of the ibidi labware, such as the [μ-Dish 35 mm, high](#) or the [μ-Slide 8 Well high](#), can be used to culture single cells in a 3D matrix, and are ideal for high-end microscopy.



Chemotaxis and Migration Assays in 3D

A [chemotaxis assay](#) is used to analyze directed cell migration towards a chemoattractant. Culturing cells in a 2D environment during a chemotaxis assay might not reflect the *in vivo* situation, resulting in an altered cell behavior and migration. To overcome this issue, cells can be embedded in a 3D matrix that mimics their natural environment, such as collagen, Matrigel®, or other hydrogels.

Advantages of 3D Chemotaxis Assays

- More *in vivo*-like setting for most cell types
- Highly defined environment (e.g., fibers or matrix)
- Chemotaxis assays with suspension cells possible

Limitations of 3D Chemotaxis Assays


- Gel handling: more parameters to control during the experiment
- Cells might attach to 2D surface, thus creating 2.5D conditions
- Cells might go out of focus during 3D tracking

Find more information about 2D and 3D chemotaxis assays in the following Application Notes:


- [AN 17: Chemotaxis 2D and 3D \(PDF\)](#)
- [AN 23: 3D Chemotaxis Protocol for Non-Adherent Cells in a Gel Matrix \(PDF\)](#)
- [AN 24: Chemotaxis of HT-1080 Cells in 2D and 3D \(PDF\)](#)
- [AN 26: Collagen I Gel for 3D Cell Culture \(PDF\)](#)
- [AN 34: Chemotaxis of HUVEC Cells in 2D and 3D \(PDF\)](#)

Biswenger V, et al. Characterization of EGF-guided MDA-MB-231 cell chemotaxis *in vitro* using a physiological and highly sensitive assay system. *PLoS One*, 2018, 10.1371/journal.pone.0203040.

[Read article](#)



Microscopy and schematic of adherent HT-1080 cancer cells on a 2D surface (left), and embedded into a 3D [Collagen I](#) gel (right) in the [μ-Slide Chemotaxis](#).



Spinning disk confocal time-lapse microscopy of LifeAct TagRFP transfected HT-1080 cancer cells, which are migrating in a 3D [Collagen](#) matrix in the [μ-Slide Chemotaxis](#), 63x oil immersion.

ibidi Solutions


The [μ-Slide Chemotaxis](#) and the [sticky-Slide Chemotaxis](#) are ideally suited for both 2D and 3D experiments. Chemotactic gradients can be easily established in water-based 3D gels, such as [Collagen I](#) gels and Matrigel®, because the gel structure does not hinder the formation of a soluble gradient by diffusion.



3D Cell Culture Under Flow


Interstitial Flow

In vivo, many cell types are constantly exposed to liquid flow. When culturing them in an *in vitro* 3D matrix, a soft interstitial flow can be applied by perfusing them with growth medium or any reagent or drug of choice. By doing this, conditions close to the cells' natural environment can be established.



Perfusion

The combination of cells inside a 3D matrix and a channel above allows for an easy application of flow. This experimental setup passively feeds the cells inside the *in vitro* 3D matrix by diffusion through the gel. Oxygen and nutrients are supplied by the gentle flow. The adjustable flow rate defines the level of nutrition, enabling long-term live cell experiments.



ibidi Solutions

The ibidi [Channel Slides](#), including the [μ-Slide III 3D Perfusion](#), the [μ-Slide I Luer 3D](#) and the [μ-Slide VI](#) families, allow for the seeding of the cells in a 3D matrix and the application of flow (e.g., using the [ibidi Pump System](#)).




Single Cells in a 3D Matrix

3D Live Cell Imaging of Migrating HT-1080 Cancer Cells in a Collagen Matrix

LifeAct-expressing HT-1080 cells (green) were seeded in a 1.5 mg/ml [Collagen Type I, Rat Tail](#) layer (white) in the [μ-Slide Chemotaxis](#). Cell migration was documented by taking a photo every 300 seconds on a Zeiss Confocal Microscope LSM 880 AxioObserver using a water immersion objective lens 40x/1.2.

[Click here](#) to watch the movie on our website.




Spheroid and Organoid Culture

Spheroid Formation on a Defined Micropattern


Micropatterns are powerful tools for optimizing 3D assays. ibidi's [μ-Slides With Multi-Cell μ-Patterns](#) enable spatially defined cell adhesion for various 2D and 3D cell culture applications. Defined adhesion spots, surrounded by [Bioinert](#), are able to catch all adherent single cells from a cell suspension. Bioinert is fully non-cell-attachable. This forces all cells to aggregate to each other at the adhesion spots, thus forming spheroids in a defined and controllable way.

Spheroid formation of NIH-3T3 cells (murine embryo fibroblasts) on 200 μm adhesion spots, Spheroid generation was documented for 64 hours. Phase contrast live cell imaging, 4x objective lens.

[Click here](#) to watch the movie on our website.




Spheroid formation of NIH-3T3 cells (murine embryo fibroblasts). Image taken 14 days after seeding single cells in the [μ-Slide 8 Well With Multi-Cell μ-Pattern](#). Brightfield microscopy, 10x objective lens.



Spheroid Formation on a Defined Micropattern


Spheroid formation of NIH-3T3 cells, 14 days after seeding single cells in the [μ-Slide VI⁹⁶⁴ With Multi-Cell μ-Pattern](#). Phase contrast microscopy, 4x objective lens.



Invasion of HT-1080 Cancer Cells in a 3D Collagen Gel

Invasive human fibrosarcoma cancer spheroids (HT-1080) were embedded into [Collagen Type I, Rat Tail](#) gel. The invasion into the gel matrix was recorded for 48 hours in the [μ-Slide 8 Well](#). 4x objective lens, brightfield.


[Click here](#) to watch the movie on our website.



Sprouting of Endothelial Cells in a 3D Collagen Gel

Live cell imaging of a spheroid of Human Umbilical Vein Endothelial Cells (HUVEC), embedded into a 3D gel made of [Collagen Type I, Rat Tail](#). The sprouting process into the gel matrix was recorded for 44 hours in the [μ-Slide 8 Well](#). 10x and 4x objective lens, brightfield.

[Click here](#) to watch the movie on our website.




Highly Improved Spheroid Growth Rates When Cultured Under Perfusion

L929 fibroblasts show spheroid formation in the [μ-Slide Spheroid Perfusion](#), Bioinert, days 1–14, seeding concentration 5×10^5 single cells/ml.


Left: no perfusion, medium exchange every second day. Right: perfusion with the [ibidi Pump System](#), 0.75 ml/min. Phase contrast microscopy, 10x objective lens, well diameter 800 μm.

[Click here](#) to watch the movie on our website.



Live/Dead Staining of a Long-Term Cultured Spheroid

Live/dead FDA/PI staining of an L929 spheroid in the [μ-Slide Spheroid Perfusion](#), Bioinert, after 14 days in culture with perfusion using the [ibidi Pump System](#), 0.75 ml/min. Green: living cells (fluorescein diacetate, FDA); red: dead cells (propidium iodide, PI). Widefield fluorescence microscopy, 10x objective lens.





Chemotaxis and Migration Assays in 3D

Chemotaxis of Endothelial Cells in a 3D Collagen I Gel Towards an FCS Gradient

Live cell imaging of Human Umbilical Vein Endothelial Cells (HUVEC) embedded in a 1.5 μg/ml Collagen Type I gel in the [μ-Slide Chemotaxis](#), migrating towards fetal calf serum. Note: the cells connecting to each other form strings during the chemotaxis process. Phase contrast, 4x objective lens, 24 hours.

[Click here](#) to watch the movie on our website.





Manufacturer

ibidi GmbH
Lochhamer Schlag 11
82166 Gräfelfing
Germany

Toll free within Germany:

Phone: 0800/00 11 11 28
Fax: 0800/00 11 11 29

International calls:

Phone: +49 89/520 46 17-0
Fax: +49 89/520 46 17-59

E-Mail: info@ibidi.com
ibidi.com

North American Headquarters

ibidi USA, Inc.
2920 Marketplace Drive, Suite 102
Fitchburg, WI 53719
USA

Toll free within the US:

Phone: +1 844 276 6363

International calls:

Phone: +1 608 441 8181
Fax: +1 608 441 8383

E-Mail: ibidiusa@ibidi.com
ibidi.com

All ibidi products are for research use only! Errors and omissions excepted.

© ibidi GmbH
FL_AG_040, V 2.1 2021/09

For free samples, application notes, and handling movies,
please visit us at:

ibidi.com