

Microfluidic slide for 3D-bioprinting applications using the example of a proximal tubule model

Philipp Linke^{1,*}, Isabelle Schmidt², Elin Pernevik³, Volodymyr Kuzmenko³, Hanna Hartmann², Zeno Guttenberg¹

¹ ibidi GmbH, Gräfelfing, Germany

² NMI Natural and Medical Sciences Institute at the University of Tübingen, Reutlingen, Germany

³ CELLINK Bioprinting AB, Gothenburg, Sweden

Introduction


Challenge: Bioprinting has come into focus as a method to create specialized organ-on-a-chip models, providing high flexibility and reproducibility of the models. However, realizing such models is usually time-consuming, costly and requires a high degree of expertise.

Our solution: A microfluidic device which makes 3D-printing on a chip accessible and allows perfusion of arbitrary hollow structures.


Methods: In collaboration with CELLINK and the NMI, we tested several

bioink candidates to develop a proximal tubule model. The rheological properties of all bioinks were assessed. Adhesion and monolayer formation of the cells were tested. The viability of the cells was assessed, and tight junction formation was confirmed by staining. We also measured the diffusion of fluorescent dextran through the layer.


Results: The cells formed a functional monolayer and created a viable model of the proximal tubule.




How to create hollow tubes




Adhesion and monolayer formation on 2D hydrogel




Viscosity and elasticity of the bioinks



Immunostaining inside tubes



Diffusion through cell layer



Summary

Results: GelMA 10 % was chosen to be used as a surrounding ink, since it facilitates optimal adhesion and monolayer formation of the cells. The printability of Pluronic as sacrificial material was characterized and fine-tuned. After creating the channels and seeding cells, their viability was shown to be sustained, and tight junction formation was clearly visible. Finally, perfusion with a fluorescent dye showed the effective barrier function of the cell layer, indicating that the cells formed a viable model of the proximal tubule.

Outlook: The universal approach of 3D-printing allows for printing of arbitrary structures, which can be adapted to mimic specific geometries of tissue models, opening a large field of applications. We believe that the benefits of our platform will enable scientist to develop tissue models more quickly and perform pharmaceutical studies on a large scale.

Acknowledgments

This work is funded by the EUROSTARS grant project PTubeChip E! 113681 and by the Federal Ministry of Education and Research (BMBF). This work also received financial support from the State Ministry of Baden-Wuerttemberg for Economic Affairs, Labour and Tourism. The authors are grateful to Martin Stelzle for the original PTube draft and Ellena Fuhrmann and Verena Singer for their technical support.