

ibidi μ -Plates

Where Precision Meets Performance

✓ Optimized for Microscopy

Brilliant optical quality and excellent cell adhesion

✓ Always in Focus

Ultra-high inner well and whole plate flatness

✓ Scalable for Every Workflow

Ranging from large growth areas to high-throughput automation (ANSI/SLAS)








Get your FREE μ -Plate at [ibidi.com/freesample](https://www.ibidi.com/freesample)



Discover the ibidi μ -Plate Family



	6 Well		24 Well		96 Well Round		96 Well Square		384 Well
Bottom									
Cat.No.	80636 80631	80637	82426 82421	82427	89606 89601	89607	89626 89621	89627	88407
Volume per well	5 ml		1 ml		200 μ l		300 μ l		50 μ l
Growth area per well	9.1 cm ²		1.54 cm ²		0.3 cm ²		0.56 cm ²		0.11 cm ²

Not finding the ideal plate for your applications?
 Contact us at: customsolutions@ibidi.com

