

# Free Sample

Use the QR Code to  
Request ibidi Labware  
Solutions for:

- Live Cell Imaging
- Migration Assays
- Immunofluorescence
- Cell Culture Under Flow
- 3D Cell Culture
- Angiogenesis Assays
- Chemotaxis Assays






Dishes



Chambered Coverslips



Chamber Slides



Channel Slides



Culture Inserts



Multiwell Plates



ibidi Mounting Media

### Choose up to 3 ibidi Products Using the ibidi Free Sample Program

Do you want to test the ibidi labware with your experiments? With the ibidi Free Sample Program, you can order free samples of the ibidi Mounting Medium and most of the ibidi labware directly from the product page.

ibidi Bottom Material Specifications	<b>POLYMER COVERSLIP #1.5</b> #1.5 ibidi Polymer Coverslip	<b>GLASS COVERSLIP #1.5H</b> #1.5H ibidi Glass Coverslip
Refractive index ( $n_D$ 589 nm)	1.52	1.52
Abbe number	56	55
Thickness	180 $\mu\text{m}$ (+10/-5 $\mu\text{m}$ )	170 $\mu\text{m}$ (+/-5 $\mu\text{m}$ )
Material	Polymer	D 263 M Schott high precision glass
Autofluorescence	Low	Low
Transmission	Very high (even ultraviolet)	High (ultraviolet restrictions)
Surface modifications	ibiTreat (tissue culture-treated) Uncoated (hydrophobic) Bioinert (surface passivation) Protein coatings	Uncoated glass