



MIDI 40

Cell Culture CO₂ Incubator
for ibidi Assays

- ✓ 40 liter volume
- ✓ Ideal for separate ECIS or perfusion experiments
- ✓ For use with the ibidi Pump System with up to 4 Fluidic Units (Fluidic Unit Quad)

Technical Features:

Outer dimensions		CO ₂ range	0–20%
Height	46.5 cm	CO ₂ pressure	1 bar (14.5 psi)
Width	47.0 cm	Voltage	250 V / 120 V (on request)
Depth	59.7 cm	Humidity	Water tray, up to 95%
Inner dimensions		Weight	27 kg
Height	35.5 cm	Stackable	Yes
Width	30.5 cm	Portable	No
Depth	35.5 cm	Access hole	25 mm
Temperature range	5°C over RT up to 60°C	Material inside	Stainless steel
Temperature uniformity	+/- 0.4°C	Door	Outer door + Inner glass door

Ordering Information:

Cat. No.	Description
10992	MIDI 40 Cell Culture CO ₂ Incubator for ibidi Assays