





FREE SAMPLE!

Go to [ibidi.com](https://www.ibidi.com) to get started!

	 POLYMER COVERSIP #1.5 #1.5 ibidi Polymer Coverslip	 GLASS COVERSIP #1.5H #1.5H ibidi Glass Coverslip
Refractive index (n_D 589 nm)	1.52	1.52
Abbe number	56	55
Thickness	180 µm (+10/-5 µm)	170 µm (+/-5 µm)
Material	Polymer	D 263 M Schott high precision glass
Autofluorescence	Low	Low
Transmission	Very high (even ultraviolet)	High (ultraviolet restrictions)
Surface modifications	ibiTreat (tissue culture-treated) Uncoated (hydrophobic) Bioinert (surface passivation) Protein coatings	Uncoated glass